

HEALTH EXPERT

Beryl®

Ultrasonic Height



- Ultrasonic to measure the height
- Measure the weight and height



We concern your health every day!

Beryl®



Model:	BYH09
Height measure way:	Ultrasonic
Weight Range:	5~200kg
Height Range:	90~195CM
Weight Precision:	division 50g,+/-300g
Height Precision:	division 0.1 cm,+/-0.5cm
LCD	68X44cm
Power	5V,0.5-3A DC adapter

How to use

Weight/height Mode:

Four-step instruction diagram for using the scale:

- Put the machine on the flat floor. Press the Mode button to weight and Height mode.
- Stand on the scale, wait for 10 second. The machine will show your weight and height on it.
- Put the machine on the flat floor. Press the Mode button to Weight/height/body.
- Stand on the scale, wait for 10 second. The machine will show your Weight/height/body fat on it.

USE RANGE

Icons representing various locations where the scale can be used:

- Hospital
- Beauty Salon
- GYM
- Chain store
- Super Market
- Hotel
- Home
- School

Health care

